

**New Kits!
RNA Promotion!
Discount 20%**

BIOTOPICS

RNA Promotion & NEW KITS • SPRING 2018



MACHEREY-NAGEL - RNA Promotion



RNA Purification - Do it right!

Tailormade solutions for different starting materials. Superior reproducibility and reliability.

NEW! NucleoBond® RNA Soil

RNA purification from soil samples for metagenomic analysis

Anion exchange technology to optimize RNA yield and purity - suitable for metagenomic studies. Parallel preparation of RNA and DNA in one hour.

NEW! NucleoSpin® RNA Stool

Speedy isolation of total RNA from various stool specimen

Combination of mechanical homogenization and chemical lysis ensures optimized RNA yield, including small RNA. Isolated RNA does not show any PCR inhibition due to NucleoSpin® Inhibitor Removal Columns.

NEW! NucleoSpin® RNA Plus XS

RNA extraction from extra small samples down to single cells analysis

gDNA removal column and a novel lysis buffer without β -mercaptoethanol or TCEP.
Elution volume from 5 μ l.

BIOTOP JL

www.biotop.fi

Itäinen Pitkätatu 4 C
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

Tilaukset ja tiedustelut

Matti Raitio
Pia Kuusisto
Risto Ranta
City-Lab Oy

puh. 050-3611 848
puh. 040-4844 854
puh. 050-5280 789

matti.rautio@biotop.fi
pia.kuusisto@biotop.fi
risto.ranta@biotop.fi

BioCity, Biomedicum, Cultivator II, ARVO, Aapistie ja Snellmania

**New Kits!
RNA Promotion!
Discount 20%**



Posti Green

RNA PROMOTION! 20% discount!

RNA from cells and tissue

	CODE	PRICE	PROMOTION
NucleoSpin® RNA Plus (50)	740984.50	287€	230€
NucleoSpin® RNA Plus (250)	740984.250	+261€	1009€
NEW! NucleoSpin® RNA Plus XS (50)	740990.50	327€	262€
NEW! NucleoSpin® RNA Plus XS (250)	740990.250	+438€	1150€
NucleoZOL (200 ml)	740404.200	267€	214€
Nucleospin® RNA set for NucleoZOL *)	740406.50	211€	169€

MicroRNA

NucleoSpin® miRNA (50)	740971.50	323€	258€
NucleoSpin® miRNA (250)	740971.250	+416€	1133€
NucleoSpin® miRNA Plasma (50)	740981.50	281€	225€
NucleoSpin® miRNA Plasma (250)	740981.250	+218€	974€

RNA, DNA, and protein

NucleoSpin® TriPrep (50)	740966.50	390€	312€
NucleoSpin® TriPrep (250)	740966.250	+699€	1359€
NucleoSpin® RNA/Protein (50)	740933.50	357€	286€
NucleoSpin® RNA/Protein (250)	740933.250	+539€	1231€

Small and large RNA from FFPE samples

NucleoSpin® totalRNA FFPE (50)	740982.50	302€	242€
NucleoSpin® totalRNA FFPE (250)	740982.250	+314€	1051€
NucleoSpin® totalRNA FFPE XS (50)	740969.50	329€	263€
NucleoSpin® totalRNA FFPE XS (250)	740969.250	+430€	1144€

RNA from plant & fungi, stool, soil

NEW! NucleoSpin® RNA Plant and Fungi (50)	740120.50	281€	225€
NEW! NucleoSpin® RNA Plant and Fungi (250)	740120.250	+235€	988€
NEW! NucleoSpin® RNA Stool (50)	740130.50	411€	329€
NEW! NucleoBond® RNA Soil (20)	740140.20	348€	278€
NEW! DNA Set for NucleoBond® RNA Soil (20)	740141.20	69€	55€

*) **Please note!** Compatible with all ZOL- and TRI-reagents!



For more
RNA Purification kits,
please check
www.mn-net.com/bioanalysis.