



BIOTOPICS

PROMOTION! Midori Green Stains and Agarose

PROMOTION!

FastGene® Agarose for DNA/RNA Electrophoresis

Molecular biology grade agarose is easy to dissolve and enables a rapid gelling. The purity of the agarose leads to an excellent transparency and a low background. DNase/RNase free.

This results in sharp and well defined DNA and/or RNA bands. Available also in 0.5 g tablets.

ORDERING INFORMATION:

AG01	FastGene Agarose (100 g)	53€	48€
AG02	FastGene Agarose (500 g)	216€	165€
AG03	Agarose (5 g) SAMPLE		FREE
AG05	FastGene Agarose Tablets (200)	107€	98€
AG09	Midori Green Agarose Tablets (75)	180€	160€

www.nippongenetics.eu



Very low
EEO activity!

PROMO!
165€
Agarose
500g

NEW!

NIPPON
Genetics

MIDORI Green Agarose Tablets

The agarose tablet contains **Midori Green DNA/ RNA Stain, TBE powder and agarose** therefore only having to be solved in pure water. After heating and pouring the solution, the gel is ready for loading.

Tablet format – no weighing equipment required.
JUST ADD PURE WATER.

Midori stain is a **non-carcinogenic alternative.**



One tablet
= 1% 50 ml
gel

MIDORI Green Agarose Tablets (AG09)

PROMOTION! +80€ 160€ (75 tablets)

www.nippongenetics.eu



Posti Green

NIPPON GENETICS

Midori Green DNA/RNA Stains

The safe alternative to ethidium bromide


ORDERING INFORMATION	PRICE
MG06 Midori Green Direct, 1 ml	99€
MG04 Midori Green Advance, 1ml	99€

Promo!
225€
3x1ml

Gel Documentation Systems with Blue/Green LED

The broad spectrum Blue/Green LED technology improves the DNA/RNA imaging sensitivity and enables the detection of Midori Green, the next generation green dyes, as well as Ethidium Bromide.

BLUE/GREEN LED + MIDORI GREEN DIRECT
= MORE SENSITIVE than UV + EtBr



The absorbance spectrum of Ethidium Bromide (red), Midori Green Advance (green) and Midori Green Direct (light green).

NIPPON Genetics

www.nippongenetics.eu

BIOTOP

www.biotop.fi

PharmaCity
Itäinen Pitkätatu 4 C
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

NIPPON GENETICS - **MUPID PROMOTION**

Mupid[®] One Electrophoresis

The Mupid™ One is an agarose gel electrophoresis system.

All accessories are included:
power controller;
gel trays, combs...



ORDERING INFORMATION

ORDERING INFORMATION	PRICE
MU2 MUPID One Electrophoresis System 1 unit	773€
MU2 MUPID One Electrophoresis System 2 units	1159€*

***PROMOTION!** Second unit discount 50%
1 unit = 773€ 2 units = **1159€** (773€+386€)

NIPPON GENETICS - DYES FOR EASY TRACKING

FastGene DNA Markers

ORDERING INFORMATION

ORDERING INFORMATION	PRICE
MWD50 50 bp 50 µg / 500 µl	47€
MWD100 100 bp 50 µg / 500 µl	47€
MWD1 1 kb 50 µg / 500 µl	36€

www.nippongenetics.eu

Tilaukset ja tiedustelut

Matti Raitio p. 050-3611 848 matti.raitio@biotop.fi
Pia Kuusisto p. 040-4844 854 pia.kuusisto@biotop.fi
Risto Ranta p. 050-5280 789 risto.ranta@biotop.fi
City-Lab Oy BioCity, Biomedicum, Cultivator II, Arvo, Aapistie ja Snellmania