



MACHEREY-NAGEL

# RNA Purification - Do it right!

Tailormade solutions for different starting materials. Superior reproducibility and reliability.

## PROMOTION!



## NucleoSpin® RNA Plus

Ultrafast and convenient

- Lysate clearing and gDNA removal with one column only
- No time-consuming rDNase digestion
- New efficient lysis buffer, **no  $\beta$ -mercapthoethanol or TCEP**

## NucleoSpin® RNA Plus XS


- RNA extraction from extra small samples **down to single cells analysis**
- gDNA removal column and a novel lysis buffer **without  $\beta$ -mercaptho-ethanol or TCEP**
- Elution volume **from 5  $\mu$ l**

For more RNA Purification kits, please check

[www.mn-net.com/bioanalysis](http://www.mn-net.com/bioanalysis)

### APPLICATION DATA

A) Reliable RNA isolation allows single cell analysis



**NucleoSpin® RNA Plus XS enables reliable isolation of high quality RNA from extra small samples down to a single cell in outstanding quality.**

A) RNA was isolated from HeLa cells (10<sup>6</sup> to 1 cell) using NucleoSpin® RNA Plus XS. RNA was analyzed by qRT-PCR.

### ORDERING INFORMATION

		PRICE	PROMO
740984.50	NucleoSpin® RNA Plus (50 preps)	287€	229€
740984.250	NucleoSpin® RNA Plus (250 preps)	+261€	999€
740990.50	NucleoSpin® RNA Plus XS (50 preps)	327€	259€
740990.250	NucleoSpin® RNA Plus XS (250 preps)	+438€	1149€

Promotion prices valid until the end of November 2018. Vat 0%.

MACHEREY-NAGEL

BIOTOP JL



BIOTOP OY  
Itäinen Pitkätatu 4 C  
20520 TURKU  
puh. (02) 241 0099  
info@biotop.fi

### Tilaukset ja tiedustelut

Biotop Oy  
PharmaCity  
[www.biotop.fi](http://www.biotop.fi)

City-Lab Oy  
Biomedicum, BioCity, Cultivator II, Arvo, Aapistie ja Snellmania  
[www.city-lab.fi](http://www.city-lab.fi)

