


NEW! Pipette Checkit™




Verify pipette accuracy in seconds

Evaluate the accuracy of your pipettes in seconds with our convenient Pipette Checkit® cartridges. Low-cost, fast, and easy to use. Available for various size pipettes.


Verification on the spot - it doesn't get any easier!




STEP 1: FILL IT
Fill pipette tip with dye



STEP 2: DISPENSE
Dispense dye into well



STEP 3: FLIP IT
Flip up well tab



STEP 4: CHECK IT!

ORDERING INFORMATION

Pipette size	Checkit® Pro model	Measured volumes (µL)		
2 µL	CP-C02-5	2	1	0.2
5 µL	CP-C05-5	5	2.5	0.5
10 µL	CP-C10-5	10	5	1
20 µL	CP-C20-5	20	10	2
50 µL	CP-C50-5	50	25	5
100 µL	CP-C100-5	100	50	10
200 µL	CP-C200-5	200	100	20

NEXT ADVANCE
Innovative Lab Products for the Life Sciences

www.nextadvance.com



For every received order, we donate to Baltic Sea protection.

BIOTOP JL
www.biotop.fi

Itäinen Pitkäkatu 4C
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

Tilaukset ja tiedustelut

Pia Kuusisto
Matti Raitio
Risto Ranta
Lasse Matintalo
City-Lab Oy

puh. 040-4844 854
puh. 050-3611 848
puh. 050-5280 789
puh. 0400-800 626
BioCity, Biomedicum, Cultivator II, Arvo, Aapistie ja Snellmania

pia.kuusisto@biotop.fi
matti.rautio@biotop.fi
risto.ranta@biotop.fi
lasse@matintalo.com