



BIOTOPICS

NEW! FAS-Digi PRO & MIDORI Green Xtra

Still destroying your DNA with UV-light?

FAS-Digi PRO

NEW! Gel Documentation System

Broad spectrum Blue/Green LED technology doesn't damage nucleic acids and allows detection of all green and red DNA dyes. Ultra-sensitive DNA signals!

FAS-Digi PRO is composed of Blue/Green LED transilluminator, imaging chamber with 24 MPixel camera and it is controlled with an imaging software on your PC.

Now **1000€ refund*** for your old UV table when purchasing new FAS-Digi PRO!

PROMOTION!

Normal price 5950€

***Price 4950€**



NEW!

MIDORI Green Xtra

Non-carcinogenic DNA stain of the latest generation - Ultra-sensitive for LED illumination



A new highly sensitive fluorescent stain for a safe visualization of DNA and RNA in agarose gels.

Agarose gels stained with MIDORI Green Xtra have a very low background fluorescence, which makes the identification of low amounts of DNA very easy.

Genetics NIPPON **Genetics**

www.nippongenetics.eu



MUPID One, MU2



Posti Green

NIPPON GENETICS - Electrophoresis



Midori GREEN DNA/RNA Stains

The safe alternative to ethidium bromide

ORDERING INFORMATION	PRICE
MG10 NEW! Midori Green Xtra, 1 ml* 99€	255€ (3x1 ml)
MG06 Midori Green Direct, 1 ml 99€	255€ (3x1 ml)
MG04 Midori Green Advance, 1ml* 99€	255€ (3x1 ml)

* Sufficient to stain 25-50 litres of agarose.

BGLED + Midori Green
= More sensitive than UV+EtBr

Midori Green dye	Blue/Green LED	Blue LED	UV-light	In-Gel	Post-staining	Direct
Advance	✓	✓	✓	✓	✓	✓
Direct	✓	✓	✗	✗	✗	✓
Xtra	✓	✓	✗	✓	✓	✗

Compatible ✓ Compatible, but not recommended ✗ Not compatible ✗

www.nippongenetics.eu

BIOTOP

www.biotop.fi

PharmaCity
Itäinen Pitkätatu 4 C
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

Mupid® One Electrophoresis System

for agarose gels. All accessories are included: power supply, gel trays, combs...

FastGene® Agarose

Molecular biology grade agarose for DNA/RNA Electrophoresis. The purity of the agarose leads to an excellent transparency and a low background. DNase/RNase free.

FastGene® DNA Markers - Dyes for easy tracking.

ORDERING INFORMATION

MU2	MUPID One Electrophoresis System	787€
AG02	FastGene Agarose (500 g)	165€
FG-830	FastGene Agarose Gel Band Cutter	41€
MWD50	FastGene DNA Marker 50 bp , 50 µg/500 µl	48€
MWD100	FastGene DNA Marker 100 bp , 50 µg/500 µl	48€
MWD1	FastGene DNA Marker 1 kb , 50 µg/500 µl	37€

Prices valid until the end of May 2019. Vat 0%

Tilaukset ja tiedustelut

Matti Raitio p. 050-3611 848 matti.raitio@biotop.fi
Pia Kuusisto p. 040-4844 854 pia.kuusisto@biotop.fi
Risto Ranta p. 050-5280 789 risto.ranta@biotop.fi
City-Lab Oy www.city-lab.fi