

**4titude
PROMOTION!**

FrameStar® PCR Plates

4titude®

FrameStar® 2-component PCR Plates (polypropylene tubes with a rigid polycarbonate frame) reduces evaporation from PCR plates, improving your results and allowing you to reduce the volume of expensive reagents to save you money.

Specific FrameStar® versions designed for all qPCR Instruments available
- with clear, white, optical or black wells. **NOTE! 50 plates/pack.**

**Buy 4 cases of
FrameStar® Plates
and get 5!**

FrameStar® 96/384 PCR Plates - 50 plates/box		Roche LC	Biorad CFX	ABi	ABiF	Eppendorf	PRICE
4ti-0720/C	96 Non-skirted, low-profile, clear wells (50)		◆				159€
4ti-0730/C	96 Semi-skirted with upstand, clear wells (50)			◆			159€
4ti-0910/C	96 FastPlate, semi-skirted, clear wells (50)				◆		159€
4ti-0770/C	96 well Semi-Skirted PCR Plate, ABI® Style			◆			159€
4ti-0960/C	96 Skirted, low-profile, clear wells (50)		◆			◆	159€
4ti-0952	96 Semi-Skirted, Roche Style, Plus qPCR Seal (50+50)	◆					274€
4ti-0384/C	384 well Skirted PCR Plate, clear wells (50)		◆	◆	◆	◆	159€
4ti-0382	384 well Skirted PCR Plate, Roche Style, Plus qPCR Seal (50+50)	◆					274€

4titude® recognises that designated trademarks and brands are the property of their respective owners.

Made in UK under ISO 13485

96 Well Standard PCR Plates

Clear polypropylene. **50 plates/box!**

**Buy 4 cases
of Standard PCR
Plates and get 5!**



Code	Product	PRICE
4ti-0740	Skirted, low profile wells	96€
4ti-0750	Non-Skirted, standard profile wells	96€
4ti-0760	Semi-Skirted, standard profile wells	96€
4ti-0736	Semi-Skirted. ABI® Style optical wells	135€

Prices valid until the end of March 2019. Vat 0%. **PROMOTION CODE SP-0609.**

Made in UK under ISO 13485

BIOTOP JL

www.biotop.fi

Itäinen Pitkäkatu 4 C
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

Tilaukset ja tiedustelut

Matti Raitio puh. 050-3611 848 matti.raitio@biotop.fi
Pia Kuusisto puh. 040-4844 854 pia.kuusisto@biotop.fi
Risto Ranta puh. 050-5280 789 risto.ranta@biotop.fi
City-Lab Oy www.city-lab.fi