

Introducing the Newest

Bullet Blender® Homogenizer Storm Pro

- with Easy Touch Technology.

The coolest way to lyse cells, homogenize tissue . . . and more

The Bullet Blender® enables you to homogenize, disrupt, or lyse up to 24 tissue or cell culture samples at a time. Load the samples into standard polypropylene tubes, then place them in the Bullet Blender. The "bullets" in the instrument vigorously strike all of the tubes and homogenize your sample in minutes. The Bullet Blender homogenizer Gold models keep samples at 4°C!



- Consistent, high yields
- Samples stay cool
- Widest array of applications, easiest to use
- No cross contamination
- Risk free

PROMOTION!

BT24M Bullet Blender Storm Pro 4800€

Homogenizer for 24 samples in microcentrifuge tubes. More powerful motor. Includes starter kit of beads, scoop and tubes. With air cooling, CE certification and touch screen.

BT5E Bullet Blender 5E Pro 4900€

Homogenizer for 12 samples in 5mL snap-cap tubes. Includes starter kit of beads, scoop and tubes. With more powerful motor, Air Cooling, CE certification and touch screen.

For more information Biotop Oy or City-Lab Oy. Promotion valid until the end of December 2020.

For every order we donate to Baltic Sea protection.



MACHEREY-NAGEL:

RNA and DNA purification from hard to lyse samples

Efficient and fast RNA and DNA extraction

NucleoSpin® Bead Tubes improve lysis efficiency.
Compatible with common disruption devices.



Orders, technical support: