

For every order we donate to Baltic Sea protection!



MACHEREY-NAGEL

Plasmid & Clean-up PROMOTION



High purity DNA isolation. Kits for transfection-grade and endotoxin-free plasmid DNA purification & clean-up for reaction mixtures and gel extraction.



20% off
plus the MN
**VIRTUAL
DUCK***

NEW! NucleoSpin® Gel and PCR Clean-up XS

2 in 1 clean-up kit for reaction mixtures and gel extraction.
XS column allows elution in ≥ 6 uL.

NEW! NucleoSpin® Inhibitor Removal

Get pure and amplifiable DNA for your downstream application!
Removes PCR inhibitors and discoloration. Superior DNA recovery.

Request a free sample kit!



*As long as stocks last. Promotion is valid until the end of June 2021.

BIOTOP JL



PharmaCity
Itäinen Pitkätatu 4 C, 20520 TURKU
puh. (02) 241 0099
info@biotop.fi
www.biotop.fi

Orders, technical support:

Matti Raitio
Risto Ranta
City-Lab Oy

puh. 050-3611 848 matti.rautio@biotop.fi
puh. 050-5280 789 risto.ranta@biotop.fi
BioCity, Biomedicum, Biokeskus 1, Aapistie ja Snellmania

MACHEREY-NAGEL

Plasmid & Clean-up PROMOTION



Application	Scale	Typical yield	Quality Grade EU (EU=endotoxin unit)	Kit name
TRANSFECTION (non-sensitive cells)	Mini	15-30 µg	10-50 EU/µg DNA	NucleoSpin® Plasmid Transfection-grade
	96-well	5-20 µg	10-50 EU/µg DNA	NucleoSpin® 96 Plasmid Transfection-grade
	Midi	500 µg	<1 EU/µg DNA	NucleoBond® Xtra Midi
	Maxi	1000 µg	<1 EU/µg DNA	NucleoBond® Xtra Maxi
TRANSFECTION (sensitive cells)	Midi	500 µg	<0.05 EU/µg DNA	NucleoBond® Xtra Midi EF
	Maxi	1000 µg	<0.05 EU/µg DNA	NucleoBond® Xtra Maxi EF
	96-well	<50 µg	<0.05 EU/µg DNA	NucleoBond® Xtra 96 EF

Endotoxin free DNA enhances transfection efficiency!

CLEAN-UP	PCR Clean-up and gel extraction		PRICE	PROMO
NEW!	740611.50	NucleoSpin® Gel and PCR Clean-up XS (50)	+04€	83 €
NEW!	740611.250	NucleoSpin® Gel and PCR Clean-up XS (250)	467€	374 €
	740609.250	NucleoSpin® Gel and PCR Clean-up (250)	467€	374 €
	740986.20	NucleoSpin® Gel and PCR Clean-up Midi (20)	232€	186 €
	740610.20	NucleoSpin® Gel and PCR Clean-up Maxi (20)	344€	275 €
NEW!	740408.50	NucleoSpin® Inhibitor Removal (50)	+96€	157 €
MINI	For rapid small scale preparation			
	740588.250	NucleoSpin® Plasmid (250)	369€	295 €
	740727.250	NucleoSpin® Plasmid EasyPure (250)	288€	230 €
	740490.250	NucleoSpin® Plasmid Transfection Grade (250)	411€	329 €
MIDI	Transfection-grade plasmid DNA			
	740410.50	NucleoBond® Xtra Midi (50)	442€	354 €
	740410.100	NucleoBond® Xtra Midi (100)	818€	654 €
	740412.50	NucleoBond® Xtra Midi Plus (50)	523€	418 €
	740420.50	NucleoBond® Xtra Midi EF (50)	665€	532 €
	740422.50	NucleoBond® Xtra Midi Plus EF (50)	791€	633 €
MAXI	Transfection-grade plasmid DNA			
	740414.50	NucleoBond® Xtra Maxi (50)	+025€	820 €
	740414.100	NucleoBond® Xtra Maxi (100)	+893€	1 514 €
	740416.50	NucleoBond® Xtra Maxi Plus (50)	+145€	916 €
	740424.50	NucleoBond® Xtra Maxi EF (50)	+587€	1 270 €
	740426.50	NucleoBond® Xtra Maxi Plus EF (50)	+773€	1 418 €
GIGA	Isolation of plasmid DNA at Gigaprep scale			
	740548	NucleoBond® PC 10000 EF (5)	537€	430 €
	740593	NucleoBond® PC 10000 (5)	385€	308 €



20% off plus the MN VIRTUAL DUCK*

*As long as stocks last.

Promotion is valid until the end of June 2021.

BIOTOP JL



PharmaCity
Itäinen Pitkätatu 4 C, 20520 TURKU
puh. (02) 241 0099
info@biotop.fi
www.biotop.fi

Orders, technical support:

Matti Raitio puh. 050-3611 848 matti.rautio@biotop.fi
Risto Ranta puh. 050-5280 789 risto.ranta@biotop.fi
City-Lab Oy BioCity, Biomedicum, Biokeskus 1, Aapistie ja Snellmania