



BIOTOPICS

Bioline - The PCR Company

SENSITIVE & FAST

SensiFAST™ cDNA Synthesis Kit

Starting from 1 pg total RNA

Highly sensitive 1st strand synthesis with unbiased and robust efficiency leading to enhanced reproducibility and data accuracy.

ORDERING INFORMATION

SensiFAST™ cDNA Synthesis Kit

BIO-65053	50 rxn	195€
BIO-65054	250 rxn	625€



INSTRUMENT COMPATIBILITY chart
at www.biotop.fi: **TECHNICAL**

FOR SYBR® GREEN AND PROBE ASSAYS

SensiFAST™ qPCR Kits

Ultra-sensitive and specific kits for a fast, broad dynamic range and highly reproducible qPCR.

SensiFAST™ qPCR mixes are designed for all qPCR platforms and they are produced, unlike other common qPCR mixes, **under ISO 13485** which makes them suitable to use as an IVD component.

SensiFAST™ Probe is designed for TaqMan® probes.

ORDERING INFORMATION

SensiFAST™ SYBR®, 2000 rxn (4x5 ml) (20x5 ml)

BIO-92020	Hi-ROX Kit	655€	2875€
BIO-94020	Lo-ROX Kit	655€	2875€
BIO-98020	No-ROX Kit	655€	2875€

SensiFAST™ Probe, 2000 rxn

BIO-82020	Hi-ROX Kit	655€	2875€
BIO-84020	Lo-ROX Kit	655€	2875€
BIO-86020	No-ROX Kit	655€	2875€

ALSO AVAILABLE! SensiFAST™ Probe One-Step Kit
for fast RT-qPCR.

Prices valid until the end of December 2021. Vat 0%.



www.biotop.fi

For every order, we donate to Baltic Sea protection.



Bioline products are produced under **ISO 13485**
and are suitable for further manufacturing use
as an IVD component.



Posti Green

MyTaq™ HS Red Mix

Routine & robust Hot-Start PCR

A versatile PCR mix with antibody-mediated Hot-Start, specifically designed for highly-specific and robust PCR. Ready-to-use 2x Mix contains ultra-pure dNTPs and red dye, allowing direct gel load!

ORDERING INFORMATION

MyTaq™ HS Red Mix, 2x

BIO-25047	200 - 800 rxn	149€
BIO-25048	1000 - 4000 rxn	658€



PCR Enzyme Guide, Bioline.

meridian BIOSCIENCE™
LIFE DISCOVERED. LIFE DIAGNOSED.

BIOTOP JL
www.biotop.fi

Itäinen Pitkätatu 4C
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

Orders, technical support

Matti Raitio tel. 050-361 1848
Risto Ranta tel. 050-5280789
City-Lab Oy www.city-lab.fi

matti.raitio@biotop.fi
risto.ranta@biotop.fi



NEW!

BAND LOST?

SimpliFi HS Mix

Amplifying ALL of your DNA for more accurate
NGS results

A high-fidelity polymerase mix using aptamer-based hot-start technology, highly suited to amplification of DNA from crude samples, for fast, inexpensive target enrichment, NGS library amplification and cloning applications.

ORDERING INFORMATION

SimpliFi HS Mix (2x)

BIO-25060	2.5 ml	214€
BIO-25061	12,5 ml	749€

Prices valid until the end of December 2021. Vat 0%.