



# CappOrigami™ Reagent Reservoirs




## Fold, pour and pipette!

V-shape ensures full usage of all reagents.

Flick-up design saves storage space.

CappOrigami™ is a reagent reservoir for multichannel pipettes. Simply flick it up by pushing the sides together. The unique construction of CappOrigami™ will form a V-shape container, which makes pipetting easier and enables picking up the very last drops of the reagent.



1. Take a flat CappOrigami™ reservoir sheet out of the pack.
2. Push the two outer sides against the middle to flick the CappOrigami™ reservoir into shape.
3. Push the end triangles in to make the CappOrigami™ reservoir ready to use.

### ORDERING INFORMATION

CA40505      CAPPOrigami reagent reservoirs 30 mL  
(8- and 16-channel pipettes), bag w/ 50 pcs

### PRICE

50€

