



FrameStar® PCR Plates

FrameStar® patented two-component design combines a rigid polycarbonate frame and thin walled polypropylene tubes - thermal expansion and sample evaporation is reduced significantly, leading to an improved consistency in PCR results.

FrameStar® 96 Well Plate



Skirted low profile PCR Plate. Designed for

Bio-Rad CFX96™ / C1000 Touch™ / S1000™. **NOW FROM STOCK!**



155€
(178€)
50 plates

ORDERING INFORMATION
FrameStar® 96 PCR Plate (4ti-0960/C).
 Skirted, low-profile, clear wells, clear frame (50 plates).

Specific FrameStar® Plates designed for ALL qPCR INSTRUMENTS AVAILABLE

With clear, white, optical or black wells. Also with barcodes. **50 PLATES/BOX!**

FrameStar® 96/384 PCR Plates - 50 plates/box	Roche LC	Biorad CFX	ABI	ABiF	Eppendorf	PRICE
4ti-0720/C 96 Non-skirted, low-profile, clear wells (50)		◆			◆	176 €
4ti-0730/C 96 Semi-skirted with upstand, clear wells (50)			◆			181 €
4ti-0770/C 96 Semi-skirted, ABI Style, clear wells (50)			◆			178 €
4ti-0910/C 96 FastPlate, semi-skirted, clear wells (50)				◆		176 €
4ti-0960/C 96 Skirted, low-profile, clear wells (50). PROMO.		◆			◆	155 €
4ti-0952 96 Semi-Skirted, Roche Style + qPCR Seals (50+50)	◆					303 €
4ti-0951 96 Semi-Skirted Roche Style (50)	◆					200 €
4ti-0384/C 384 well Skirted, clear wells (50)		◆	◆	◆	◆	183 €

Now from stock!

Promotion is valid until the end of May 2022.

4titude® recognises that designated trademarks and brands are the property of their respective owners.