

BIOTOP

Cell Culture

PromoCell - The Cell Culture Experts

Cancer Cell Culture

PromoCell cancer cell culture portfolio provides innovative solutions for the isolation and culture of primary tumor cells as well as established cancer cell lines. The serum-free and xeno-free formulation of cancer cell culture media allow for an optimally adapted and highly standardized environment for 3D spheroids and 2D adherent monolayers.

Cancer Cell Culture
3D Tumorsphere Medium XF
Cancer Cell Line Medium XF
NCCD Reagent
Primary Cancer Culture System



Cell Culture Medium

High-quality medium is vital for successful cell culture. Therefore, PromoCell adheres to strict quality standards and produces its medium only from raw materials from strictly approved sources.



Growth Medium (Ready-to-use)

Consists of **Basal Medium** and **Supplement-Mix** (on vial). Complete medium is ready after one pipetting step.



Growth Medium Kit

Consists of **Basal Medium** (500 ml) and **SupplementPack**. Growth supplements come in separate vials. Allows customization of the end concentrations of growth supplements.

Human Primary Cell Culture

PromoCell supplies everything you need to cultivate human primary cells – from defined media optimized to enhance cell growth, to high quality human primary cells.

Human primary cells are available from a wide variety of tissues and organs.

Human Primary Cell culture
Chondrocyte Cell Culture
Endothelial Cell Culture
Epithelial Cell Culture
Fibroblast Cell Culture
Follicle Dermal Papilla Cell Culture
Keratinocyte Cell Culture
Melanocyte Cell Culture
Myocyte Cell Culture
Osteoblast Cell Culture
Pericyte Cell Culture
Preadipocyte Cell Culture



Human Stem & Blood Cell Culture

PromoCell offers a large range of adult stem cells and differentiated blood cells, as well as media that mimic the in vivo environment or stimulate cell differentiation.

Human Stem and Blood Cell
Dendritic Cell Culture
Hematopoietic Stem Cell Culture
Macrophage Cell Culture
Mesenchymal Stem Cell Culture
Monocyte Cell Culture
Mononuclear Cell Culture



Mycoplasma Detection, prevention & treatment



Disinfectant Solution for Incubator Water Tray

Bioguard-A Solution contains a disinfectant that does not cause damage to the stainless steel tray, is non-toxic, non-volatile, and extremely effective. We recommend the use of Bioguard-A in all water trays of common incubator types.

Disinfectant Spray for Incubators & Sterile Benches

Bioguard-S solution prevents contamination and growth of fungi (and spores), bacteria (and spores), mycoplasma, and viruses (including HIV and Hepatitis B) on surfaces. Bioguard-S is non-toxic and biodegradable.



For Rapid Detection of Mycoplasma

MycoKit is designed to detect the presence of mycoplasma in cultured cells with PCR.

ORDERING INFORMATION

Code

D1020-100
D1010-1000
K1000

Product

Bioguard-A
Bioguard-S
MycoKit PCR test

Size

100 ml
1000 ml
20 tests



biowest
The Serum Specialist



Curio C-Slide

Disposable Cell counting chamber / Hemocytometer.
Cost-saving 4-Channel slide. Neubauer Improved.
Use manually or with automatic cell counter.

FACSCOPE™ B

Fast and accurate Image-Based Cell counter.
Auto-focusing with advanced analytical algorithm using 4-channel disposable slide (Curio C-Slide). User changeable setup modes.

CURIOSIS



Biowest serum & BSA

Wide range of sera and media for cell culture, including bovine sera of different geographical origins.

FBS Premium

A selection of best serum batches, selected on excellent values of Endotoxin, Hemoglobin, and Cell growth.

Serum & BSA (Biowest)		
FBS Premium	500 ml	S181B-500
FBS Tetracycline Free	500 ml	S181T-500
Human Serum, AB male HIV testes	100 ml	S4190-100
Also available ANIMAL SERA: Bovine, Calf, Horse, Donkey, Rabbit, Goat, Lamb etc.		
Biowest BSA Bovine Serum Albumin	100 g	P6154-100
Biowest BSA Bovine Serum Albumin	500 g	P6154-500

biowest
The Serum Specialist